DIEMET-e

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E

Instructions for use of DIEMET-E

- Epoxy resins must be dosed very carefully.
- Hang the measuring/mixing bowl into the bowl support.

Dosage

- Press resin and hardener into the measuring/mixing bowl until the same graduation lines on the 2 dosing syringes are reached (the stoppers in the syringes are separated from the slides, in order to avoid influx of air if the slides are drawn back). Take drops adhering to the dosing syringes off into the bowl. Each graduation line of resin and hardener corresponds to the quantity necessary to obtain app. 4 dies, with the quantity of 3 graduation lines about one full dentition can be obtained.
- Add as many spoons full of filler to the liquid as scale graduations have been exceeded by the liquid's level on the measuring bowl. Note: If the liquid level lies between 2 graduation marks you may round up or round down. If rounded up, the viscosity of the mix augments slightly, if you round down, the resin is slightly more fluid. However, the physical properties of DIEMET-E will not alter in either case.
- Spatulate well; afterwards pour the mix, in a thin jet, into the impression.
- Clean the measuring/mixing bowl well after hardening process.

Technical data

Working time	15 min.
Gel time	20 - 30 min.
Hardening after	6 - 8 hours*
Final hardness	> 200 N/mm ²
Contraction	0,05 % or 5 µ/cm
* 90 % final hardness after	5 hours at 40° C 6 hours at 30° C 8 hours at 20° C

Hardening should be carried out at a temperature of 20° - 40° C.

Important notice:

Condensation and addition cured silicone and polysulfide impressions do not have to be insulated; remove eventual moisture from the impression and cast (addition cured silicone materials can be treated with EKM powder ord. no. 758 000 or other separating media offered by the manufacturers in order to avoid the formation of air bubbles during casting).

Polyether impressions must be insulated. This is done best with EKM powder, thus an insulating layer is avoided. The impression is dried and covered with EKM powder with a soft brush and then blown out and cast (a scarcely visible layer of powder remains). EKM powder is surrounded by DIEMET-E, it insulates reliably, but does not produce a thick insulating layer.

Hydrocolloid impressions are dried carefully (no strong air jet!) and cast immediately with DIEMET-E. Ensure that the hydrocolloid does not dry out during hardening or put it into a pressure flask for 1 hour at 1.5 bar pressure.

Seal the impressions against flowing out, as the material is very fluent.

Resin and hardener: Xn are a health danger if inhaled, swallowed and on skin contact (sensitization possible).
Do not empty these products into the soil or sewerage. Wear appropriate protective clothing.
Resin contains: 1,4 Bis (2,3 epoxypropoxy) butan
Hardener contains: benzyl alcohol phenol

 ERKODENT®

PRECISION AND PROGRESS