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## Hints

- Wounds on the feet should **always** be covered by a wound cover foil.
- In order to obtain a good level mould, a thin layer of ALGIFORM can first be filled into the casting tray (app. 500 ml).
- When using side trays the base tray should always be filled up to the brim.
- The impression should always be made under partial load.
- Protect the bottom of the tray against LEISTON with a rubber sheet.
- Use 50 g more ALGIFORM per litre of water for impressions that are filled lateron to avoid shrinkage of the negative.
- The water temperature should never exceed 25 °C.

### IMPEDO-technique

Rational moulding method for the production of negatives for inner soles and resin lasts. Skin tolerated impression material, particularly suited for patients suffering from diabetes and rheumatism.

IMPEDO trays are supplied in two lengths (30 / 35 cm) and two heights (6 / 11 cm). A heel height of 1 cm is given. ALGIFORM is supplied in packing units of 5 kg and LEISTON last resin in 5 ltr containers.

1. The foot is placed in an IMPEDO tray of suitable size. The heel height can be changed by inserting a wedge into the tray. Fill water in the base tray to find out the required capacity. 2 measuring beakers ALGIFORM (app. 100g) are mixed with 1 ltr. of water.

2. The mixed ALGIFORM is poured into the IMPEDO tray and the patients foot is placed in the base tray. In doing so it should be ensured that the leg position is at 90° and the impression should always be made under partial load. The base tray should always be filled up to the brim.

3. After ALGIFORM has hardened (app. 3-4 min.) the side trays are fit and closed with closing buttons. The second ALGIFORM mixture is prepared and cast up to the required height of the upper.

4. After hardening of ALGIFORM, the closing buttons are removed and ALGIFORM cut through at the front and rear with a separating spatula. Remove the side trays and the foot from the base tray. Check the impression for any broken out parts. The trays are rejoined and closed with the closing buttons.

5. Before pouring with LEISTON, any accumulated water is removed from the ALGIFORM impression. The required quantity of LEISTON is mixed, pour out the ALGIFORM impression and allow to harden (app. 30 min.)

6. The finished LEISTON last is removed from the ALGIFORM impression and the material further finished. LEISTON can be milled, ground, nailed, trowelled, glued and with commercial glues additions can be put on the LEISTON.

## Lasts with REVO - Gel

7. For the production of uppers for beds and arch supports made of REVO-Gel, the fore foot region is cut out to the appropriate size in the form of an arch support.

8. REVO-Gel is filled into the melting device and melt at app. 165 °C (position 10 control knob). Stir up REVO-Gel from time to time. Since heating up of the REVO-Gel takes some time, it is advisable to switch on the melting device well in advance, in order to avoid long waiting times.

9. After completely filling the negative with REVO-Gel it should be allowed to cool down thoroughly before the blank is removed from the ALGIFORM impression.

10. The finished REVO-Gel blank can be ground with commercial grinding machines. Parts which have been ground or cut off can be put back in the melting device. After some time the colour of REVO-Gel changes and with growing soilings also its properties.